

## IMPORTANT NOTICE

### REGARDING LUCIRA COVID-19 TEST KIT EXPIRY DATE


Thank you for your order. We want to provide you updated information about the expiration date on your Lucira COVID-19 Test Kits.


**The shelf life for Lucira COVID-19 Test Kits has been extended to 18 months.** For COVID-19 tests, manufacturers are required to label products based on real-time stability studies. As of June 15th, 2022, Lucira Health completed the real-time stability testing to support an 18-month shelf life.


**Lucira COVID-19 Test Kits boxes shipped from Lucira after 29 August are labeled to reflect the 18 month shelf-life. \***

**To determine the expiration date for Lucira COVID-19 Test Kits received prior to 29 Aug 2022:**

1. Locate the label below on the kit box.



2. Add 6 months to the date next to the  symbol.

The date next to  is expressed as YYYY-MM-DD where YYYY is the year, MM is the month and DD is the day.

Example: 2023-02-23 (23 February 2023) - New Expiration Date (+ 6 months): 23 August 2023.

\*Pouched items enclosed within the test kit box may have an expiration date earlier than the updated kit expiry date. However, the expiration date of these pouched items now coincides with the expiry date for the kit.

If you have any questions, please call us at 1-888-582-4724 (toll free for United States and Canada), email us at [CS@lucirahealth.com](mailto:CS@lucirahealth.com) or via our website at [lucirahealth.com](http://lucirahealth.com).

Thank you,

The Lucira Team